

MVAP

MEDICAL SUPPLIES INC

Speak to a Live Representative at:
1-877-735-6827 mvapmed.com



ORDER YOUR SUPPLIES ONLINE
MVAPMED.COM

ZBrand® Teardrop Cloth Electrode

Vermed ZBrand® Electrode line was developed to meet the requirements of the entire hospital by offering an economic, generic patient monitoring solution with quality performance.

Pack of 50

NEW!



Z05MC-50 40x36mm, 50 per pack..... \$5.15

ZBrand® Teardrop Foam Electrode

Vermed ZBrand® Electrode line was developed to meet the requirements of the entire hospital by offering an economic, generic patient monitoring solution with quality performance.

Pack of 30

NEW!



Z05SG-30 40x36mm, 30 per pack..... \$2.95

ZBrand® Teardrop Tape Electrode

Vermed ZBrand® Electrode line was developed to meet the requirements of the entire hospital by offering an economic, generic patient monitoring solution with quality performance.

Pack of 50

NEW!



Z05MT-50 40x36mm, 50 per pack..... \$5.15

ZBrand® Rectangle Foam Electrode

Vermed ZBrand® Electrode line was developed to meet the requirements of the entire hospital by offering an economic, generic patient monitoring solution with quality performance.

Pack of 50


NEW!



Z35-50 40x35mm, 50 per pack..... \$4.75

ICECap™

The ICECap™ Net uses next-generation technology to provide an easy-to-apply, disposable solution for EEG diagnostic testing anywhere, anytime. Available in sizes for adult and pediatric patients, the full 10/20 application is lightweight, flexible and adaptable to all current EEG devices. The construction is durable and lasts up to 72 hours.



NEW!

1001-01014	Adult	\$95.00
1001-01015	Pediatric.....	\$95.00
1001-15024	Adapter	\$285.00

Multi-Purpose Wraps

Velcro wraps great for securing oximeter probes or organizing EEG leads during lab procedures. 5.25" x 1.25".

Pack of 12

NEW!



MPW-12 pack of 12 \$8.95

SleepSense® Pressure Sensor/Alice 6 Compatible

1 Year Warranty



SLP14860 D/C Pressure Sensor, Alice 6 Compatible..... \$395.00

SleepRes V-Com®

V-Com is a simple product that provides comfort on inhalation for patients using PAP therapy. The high flows and pressures of CPAP that are difficult to tolerate are softened with V-Com. It's designed and engineered to soften inspiratory flow and pressure while preserving therapy (EPAP). V-Com simply makes CPAP more comfortable, which is great for new patients and those that have been on CPAP for years!



NEW!



CAP1026 SleepRes V-Com..... \$25.00

SleepSense® NightOne Compatible Reusable Belt

1 yr Warranty



SLP9010N1 \$220.00

Wire Attached Electrodes

1" round electrode with attached lead. 1.5mm connector. Sold by the pack of 3.

Packs of 3

NEW!



G60AS	60" Leads	\$5.25
G96AS	96" Leads	\$6.75
G120AS	120" Leads	\$7.50

SleepSense® Alice 6 Compatible Pediatric Thermistor


6 Month Warranty



SLP1473 \$165.00

SleepSense® PVDF Effort Kit / Alice 6 Compatible

Using PVDF technology, the Respiration Effort Sensor converts chest or abdominal respiration movement to a small analog voltage that provides a clear, reliable indication of respiration waveforms.



SLP1394	PVDF Sensor Only, Alice 6/Keyhole	\$270.00
SLP1394-KIT	PVDF Effort Kit (Includes Belt Strap), Alice 6/Keyhole	\$295.00

SleepSense® Disposable Inductive Kits

- SleepSense Disposable Inductive Belts are pre-cut into Adult sizes
- Reusable Inductive Interface Buckles have a 1 year warranty
- Easy to use, full adjustment at any time throughout the study



Kits Include: Reusable Buckles and 20 Pack of Disposable Belts

Adult Kits

SLP9107-KIT	Safety DIN Connectors, Adult 120cm	\$430.00
SLP9107A6-KIT	Alice 6 Compatible, Adult 120cm	\$430.00
SLP9207-KIT	Key Connector, Adult 120cm	\$430.00
SLP9107XL-KIT	Safety DIN Connectors, Adult Large 150cm	\$435.00

Pediatric Kits

SLP9107P-KIT	Safety DIN Connectors, Pediatric 60cm	\$430.00
SLP9154A6P-KIT	Alice 6 Compatible, Pediatric 60cm	\$430.00
SLP9207P-KIT	Key Connector, Pediatric 60cm	\$430.00

HST Kits

SLP9256-KIT	Alice NightOne Compatible, Adult 120cm	\$220.00
SLP9255-KIT	NOX T3 Compatible, Adult 120cm	\$380.00
SLP9107NO-KIT	Nomad, Adult 120cm	\$445.00

SleepSense® Alice 6 Compatible Body Position Sensor

These DC body position sensors were designed with Vertical Velcro® tabs so you can apply it directly on



SLP1584	\$270.00
---------	-------	----------

SleepSense® Alice 6 Compatible Snore Sensor

Provide clear, artifact-free qualitative indication of snoring and respiratory sounds. Wire length = 82 inches

6 Month Warranty



SLP1267	\$145.00
---------	-------	----------

EEG-acp Adhesive Conductive Paste

Adhesive Conductive Paste for Routine EEG, Sleep, Long-Term EEG and Epilepsy EEG.

- FDA Approved
- Skin Friendly
- Short Drying Time
- Excellent and Stable Impedances
- Economical in the Installation

Pack of 10



CR-050-10	Pack of 10, 100g Tubes	\$114.95
-----------	------------------------------	----------

SafePrep

- Single-use, dye free degreasing gel for EEG skin preparation
- Minimizes the risk of cross-contamination
- Proven high performance in reducing skin impedance
- Dye-free formula to minimize irritation and reduce risk of allergy in sensitive skin
- Innovative recyclable, plastic-free primary and secondary packaging



NPREP010E-1	10 gram sticks, 100 per box	\$98.00
-------------	-----------------------------------	---------

Telemetry Pouches

- Measures 6" x 7.5"
- 2.5" slit for lead wires of transmitter
- Re-closable adhesive tab secures top flap
- 5" x 64" polyester straps that can easily be tied around neck or waste

NEW!



00905A	\$2.50
--------	-------	--------

CleanCaseRx Disposable/Single Patient Use HST Case

Follow AASM HST COVID-19 guidance for using 'single use, fully disposable' testing components.

- Reduce the risk of cross-contamination between HST studies
- Build peace-of-mind in your patients and staff with one sanitary solution
- Protect the health and safety of your medical staff, patients, and HST devices
- Eliminate any need to clean, disinfect or replace expensive reusable cases
- Preserve valuable staff time and expense resources
- Accommodates FOUR HST devices and accessories in ONE CleanCaseRX - NOX Medical (T3, T3x, A1) and Itamar WatchPat
- Eco-Friendly and 100% Recyclable
- Quick and easy assembly
- Space saving and collapsable design for convenient storage



Case of 50

CCR-101 CleanCaseRx, Case of 50\$289.50

10/20 Brainet™ Templates

Individually packaged color-coded electrode positioning template with openings for 19 recording electrodes, plus openings for a ground and a reference electrode. Designed to aid in rapid and accurate placement of the full 10/20 system of scalp electrode placement.

Latex Free



BN-106	Extra Small Template (45.5-49cm)	\$63.00
BN-107	Small Template (49.5-54cm)	\$63.00
BN-108	Medium Template (54.5-57cm)	\$63.00
BN-109	Large Template (57.5-62cm)	\$63.00

Baby Brainet™ Templates

Individually packaged color-coded electrode positioning template with openings for 13 recording electrodes, plus openings for a ground and a reference electrode. Designed to aid in rapid, accurate EEG set-up in the neonatal ICU.

Latex Free



BN-106	Extra Small Template (45.5-49cm)	\$63.00
BN-107	Small Template (49.5-54cm)	\$63.00
BN-108	Medium Template (54.5-57cm)	\$63.00
BN-109	Large Template (57.5-62cm)	\$63.00

RespiPack From Dymedix

The RespiPack includes state-of-the-art, PVDF-based, reusable respiratory sensors that provide increased sensitivity and response time. The sensors are easy to clean and are compliant with AASM technical specifications.



NEW!

- Each Adult Kit includes: Reusable Airflow Sensor, Reusable Snore Sensor, (2) PerfectFit 2 or PerfectFit 2 R/S Sensors, (2) 45" Belt Straps, (2) 60" Belt Straps, FM2 Airflow Interface Cable, AccuSnore Interface Cable.

RP-001-K	Adult Kit, Alice 6	\$900.00
RP-002-K	Adult Kit, Natus Key Connector	\$900.00
RP-003-K	Adult Kit, Graef	\$900.00
RP-004-K	Adult Kit, Universal 1.5mm	\$900.00

- Each Pediatric Kit includes: Reusable Airflow Sensor, Reusable Snore Sensor, (2) PerfectFit 2 or PerfectFit 2 R/S Sensors, (2) 15" Belt Straps, (2) 25" Belt Straps, (2) 35" Belt Straps, FM2 Airflow Interface Cable, AccuSnore Interface Cable.

RP-005-K	Pediatric Kit, Alice 6	\$900.00
RP-006-K	Pediatric Kit, Natus Key Connector	\$900.00
RP-007-K	Pediatric Kit, Graef	\$900.00
RP-008-K	Pediatric Kit, Universal 1.5mm	\$900.00

TEN20™ Conductive Paste Single-Use Cups

by Weaver & Co.

Box of 24 Single Use Cups. 15g per cup.



10-20-1	1-7 Boxes	\$27.50
	8-23 Boxes	\$27.10
	24+ Boxes	\$26.70